DESCRIÇÃO

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Osteosynt®

HA/TCP micromacro porous biphasic calcium phosphate Bioceramic

Bioceramics for bone reconstruction
NCM: 30064020
Register nº Anvisa : 10273030001
(Brazilian Ministry of Health)

DESCRIPTION

OSTEOSYNT® is a third generation bioceramics composed mainly by HA (hidroxiapatite) and (βTCP) calcium (tri) phosphate presenting two phases - one more soluble and the other more stable - which define it as being biphasic and containing micro and macro pores. These characteristics give this bioceramics its properties of intrinsic osteoinduction (chemotaxis) and osteoconduction capacity (haptotaxis) required for the process of bone reconstruction. Its components and characteristics are found in both the bone mineral matrix and the dental enamel.

OSTEOSYNT® is fully biocompatible and shows the mechanical resistance necessary and desirable for the support and scaffolding of the bone reconstruction process, as it combines with the new tissue and is gradually replaced by the continuous natural process of bone remodelling.

Scientific research and clinical observations have shown the absence of rejection, undesired inflammatory, cytotoxic, immunogenic and systemic reactions. Due to its high level of calcium (Ca), which facilitates its identification, **OSTEOSYNT®** is radiopaque.

OSTEOSYNT® forms the desirable bone support matrix during the length of time required for the deposition and maintenance of the bone tissue, which firmly adheres to the surface of the matrix through a chemical process and because of the matrix surface characteristics. The bone tissue then enters its porous structure, which fully combines with the new tissue, becoming an integral part of it. **OSTEOSYNT®** acts fast and with the necessary mechanical resistance as an intrinsic conductor and inducer of the osteoblastic mitosis process of the cell by liberating calcium ions.

INDICATIONS

OSTEOSYNT® is an elective bioceramics intended for the filling and reconstruction of bone defects, as well as for the maintenance of anatomic structures, in general being indicated for use in the fields of orthopaedics and traumatology, neurosurgery, mouth-jaw-facial surgery, reconstructive surgery, craniofacial surgery, ophthalmology, bone loss and injury or anatomic structure reconstruction of segmental bone losses, deepening and dehiscence; pseudoarthrosis, bone infectious processes, osteomyelitis, osteolysis, cysts and tumors; laminectomy, arthrodesis;

PRECAUTIONS

Surgical use of **OSTEOSYNT®** must be restricted to qualified professionals with expertise in the surgical techniques and procedures indicated once the inadequate application of the product may not reach the results expected and/or trigger migration of the product. Presurgical evaluations, adequate selection of materials, use of compatible surgical techniques and procedures as well as post-surgical control are crucial to produce the desired results.

STERILE ETHYLENE OXIDE (ETO). SINGLE USE. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT USE IF CONTAINER SEAL IS BROKEN OR AFTER EXPIRATION DATE.

Unused products opened during surgery must be immediately disposed of in accordance with the current legislation about hospital waste. Out-of-date and unused products with a broken seal must be destroyed before disposal.

STORE AT ROOM TEMPERATURE

WARNING

IMPORTANT: three conditions must be present to ensure new bone formation:

- 1 - Scaffolding;
- 2 - Morphogenetic proteins (BMP's);
- 3 - Specific cells (mother cells).

Absence of any of those conditions shall compromise the results.

Complications such as particle extrusion, migration of granules, dehiscence of soft tissues, delayed scarring associated to failure to apply surgical techniques and procedures do not necessarily compromise clinical results.

This product is radiopaque, visible under X-ray.

To patients: Physical activities are restricted during the recovery period and must be evaluated according to the type of surgery, extent of lesion and area of application of the product. Deambulation may take from three to six months, provided that professional prescriptions concerning the length of time required for natural organic neo-formation of the tissue are followed.

COUNTER-INDICATIONS

OSTEOSYNT® is counter-indicated in the presence of non-removed infected and/or necrotic tissues. Its use in patients with systemic illnesses such as diabetes mellitus, AIDS, illnesses or conditions which lead to bone demineralization, or who are under corticotherapy or radiotherapy treatment does not imply undesirable reactions. **OSTEOSYNT®** does not induce immunoallergic or cytotoxic reactions. However, under those circumstances, and due to systemic implications specific of the patient, the results may not be the ones expected.

PROCEDURES

OSTEOSYNT® must be prepared in site, applied to the healthy bleeding bone after cleaning and removal of the defective bone. It should preferably be combined with the patient's autogenous medullar blood and/or autogenous medullar bone. The product can also be used pure or isolate. **OSTEOSYNT®** can also be used in combination with fibrin, platelet concentrates, cell concentrates and blood products, and/or adequate binding substances which allow its proper modelling according to indication. For its characteristics and physicochemical properties the product can be used as a conductor and releaser of such drugs as antibiotics, proteins, chemotherapeutics and others.

OSTEOSYNT® characteristic is to form a BONE MORPHOGENETIC COMPLEX.

For a faster recovery and better results the following is recommended:

Association of medullar blood puncture (up to 3.5 cm³ or 3.5ml each puncture) with **OSTEOSYNT®** in the 5g bottle. Shake the bottle for some four or five minutes. One gram of **OSTEOSYNT®** immediately absorbs approximately 0.7cm³ of medullar blood, which contains all the autogenous ingredients necessary for the process to complete. Preferably puncture the iliac or femur as they contain a higher level of the substances necessary for the process.

The process requires a compact filling of the area, as well as its adequate containment and stabilization.

PRESENTATION

Particles, parts, 10 µ micrometers micropores and macropores up to 400 µ micrometers, which characterize the product as micromacro porous.

Highly pure material presented in the following forms: granular with various diameters; blocks and parts of various forms; plastic, mouldings and coatings for use according to specifications. The product is also found in specific prototyping and 3D modelling ceramic parts or cementing composite moulding.

Cervical Block (mm)		Wedges (mm)		Bioceramic Covering with Pins for External Fixation		Bucco-maxilar Blocks (mm)	
06x04x10x12	OSTC 6.4.10.12	6x8x10x12	OSBC 6.8.10.12	200mm-M6	OSPF M6 L200	08x15x25	OSBM 8.15.25
06x04x12x12	OSTC 6.4.12.12	6x8x12x20	OSBC 6.8.12.20	200mm-M4,5	OSPF M4,5 L200	10x10x10	OSBM 10.10.10
06x04x14x12	OSTC 6.4.14.12	8x10x15x25	OSBC 8.10.15.25	150mm-M3,5	OSPF M3,5 L150	10x10x40	OSBM 10.10.40
08x06x10x12	OSTC 8.6.10.12	3x5x15x10	OSBC 3.5.15.10	100mm-M2,5	OSPF M2,5 L100	6x15x50	OSBM 6.15.50
08x06x12x12	OSTC 8.6.12.12	3x7,5x15x20	OSBC 3.7½.15.20	60mm-M2,0	OSPF M2 L60		
08x06x14x12	OSTC 8.6.14.12	3x10x15x25	OSBC 3.10.15.25				
10x8x10x12	OSTC 10.8.10.12	3x12,5x15x35	OSBC 3.12½.15.35				
10x8x12x12	OSTC 10.8.12.12	3x15x15x45	OSBC 3.15.15.45				
10x8x14x12	OSTC 10.8.14.12						

Use in Dentistry (µ)						
	10 - 20	20 - 40	40 - 60	60 - 80	80 - 100	100 - 200
0,5g	OSGD ½ [10:20]	OSGD ½ [20:40]	OSGD ½ [40:60]	OSGD ½ [60:80]	OSGD ½ [80:100]	OSGD ½ [100:200]
1,0g	OSGD 01 [10:20]	OSGD 01 [20:40]	OSGD 01 [40:60]	OSGD 01 [60:80]	OSGD 01 [80:100]	OSGD 01 [100:200]

Vials (µ)			Blocks (mm)		Buttons		Porous Sphere		
	20 - 40	40 - 60	60 - 80						
02g	OSGP 02 [20:40]	OSGP 02 [40:60]	OSGP 02 [60:80]	08x10x12	OSBP 08.10.12	10	OSBT Ø10mm	14	OSEP Ø14mm
05g	OSGP 05 [20:40]	OSGP 05 [40:60]	OSGP 05 [60:80]	05x10x20	OSBP 05.10.20	12	OSBT Ø12mm	16	OSEP Ø16mm
10g	OSGP 10 [20:40]	OSGP 10 [40:60]	OSGP 10 [60:80]	08x12x20	OSBP 08.12.20	14	OSBT Ø14mm	18	OSEP Ø18mm
				10x15x25	OSBP 10.15.25	16	OSBT Ø16mm	20	OSEP Ø20mm

Designated Responsible Technician:
Francisco Henrique Lanna Wykrota • CRM-MG 7182

Eincobio Biomaterial Ltda
CGC 00.332.420/0001-75
Av. André Cavalcanti, 63 - Gutierrez
Belo Horizonte • MG • 30430-110
www.eincobio.com.br

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Cilindro 10x15 mm OSRP ©10x15